

Incubators

Reach-In Large Capacity



Optimized Floor Space

Models SMI31 and SMI39 incubators provide 30.8 & 38.6 cu.ft. chamber capacities while minimizing the amount of floor space used. Both models incorporate our microprocessor controller to achieve precise temperature uniformity. An independent, secondary temperature controller offers the added security of over temperature protection. The chamber floors are reinforced to support roller apparatus or shakers. Both models are supplied with six shelves capable of supporting the weight of benchtop instruments.

FEATURES

- Four (4) Interior Electrical Outlets
- Prewired Panel for Chart Recorder
- Optional Reversible Door
- Tempered Glass Viewing Window

Applications

- Hematological Studies
- Microbiological Determinations
- Pharmaceutical Stability Assays
- General Purpose Incubator Applications
- Large Scale Roller Apparatus Applications
- Bacterial Culturing and Research
- Food Processing Quality Control
- Biochemical Studies



